

# PRACTICAL EXAM



Making science together!

2019-07-24



MINISTÈRE  
DE L'ÉDUCATION  
NATIONALE ET  
DE LA JEUNESSE

MINISTÈRE  
DE L'ENSEIGNEMENT SUPÉRIEUR,  
DE LA RECHERCHE  
ET DE L'INNOVATION

## General instructions

- This practical booklet contains 27 pages.
- Before the start of the practical exam, the **Read** command is given. You will have 15 minutes to read the exam booklet. You may only **read** during this time; **do not write or use the calculator**.
- You may begin working as soon as the **Start** command is given. You will then have **5 hours** to complete the exam.
- You may work on the tasks in any order, but **starting with problem P1 is advised**.
- All results and answers must be clearly written **in pen in their respective designed areas** on the exam papers. **Answers written outside the answer boxes will not be graded**.
- If you need scrap paper, use the backside of the exam sheets. Remember that **nothing outside the designed areas will be graded**.
- The official English version of the exam booklet is available upon request and serves for clarification only.
- If you need to leave the laboratory (to use the restroom or have a drink or snack), raise the appropriate card. A lab assistant will come to accompany you.
- Shelves above the benches are not to be used during the task for the purpose of equality.
- You must **follow the safety rules** given in the IChO regulations. If you break the safety rules, you will receive only one warning from the lab assistant. Any safety rule violation after the first warning will result in your dismissal from the laboratory and the nullification of your practical examination.
- Chemicals and labware will be refilled or replaced without penalty only for the first incident. Each further incident will result in the deduction of 1 point from your 40 practical exam points.
- The lab supervisor will announce a 30 minutes warning before the **Stop** command.
- You must stop your work immediately when the **Stop** command is announced. Failure to stop working or writing by one minute or longer will lead to disqualification of your practical exam.
- After the **Stop** command has been given, the lab supervisor will come to sign your answer sheet.
- After both the supervisor and you sign, place this exam booklet in the envelope and submit it for grading together with your product and thin-layer chromatography (TLC) plates.

## Lab rules and safety

- You must wear a lab coat and keep it buttoned up. Footwear must completely cover the foot and the heel.
- Always wear safety glasses or prescription glasses when working in the lab. Do not wear contact lenses.
- Do not eat or drink in the lab. Chewing gum is not allowed.
- Work only in the designated area. Keep your work area and the common work areas tidy.
- No unauthorized experiments are allowed. No modification of the experiments is allowed.
- Do not pipette with your mouth. Always use a pipette filler bulb.
- Clean up spills and broken glassware immediately from both the bench and the floor.
- All waste must be properly discarded to prevent contamination or injury. Water solutions are eligible for sink disposal. Organic waste must be disposed of in the marked capped container.

## Physical constants and equations

In these tasks, we assume the activities of all aqueous species to be well approximated by their respective concentration in  $\text{mol L}^{-1}$ . To further simplify formulae and expressions, the standard concentration  $c^\circ = 1 \text{ mol L}^{-1}$  is omitted.

Avogadro's constant:

$$N_A = 6.022 \cdot 10^{23} \text{ mol}^{-1}$$

Universal gas constant:

$$R = 8.314 \text{ J mol}^{-1} \text{ K}^{-1}$$

Standard pressure:

$$p^\circ = 1 \text{ bar} = 10^5 \text{ Pa}$$

Atmospheric pressure:

$$P_{\text{atm}} = 1 \text{ atm} = 1.013 \text{ bar} = 1.013 \cdot 10^5 \text{ Pa}$$

Zero of the Celsius scale:

$$273.15 \text{ K}$$

Faraday constant:

$$F = 9.649 \cdot 10^4 \text{ C mol}^{-1}$$

Watt:

$$1 \text{ W} = 1 \text{ J s}^{-1}$$

Kilowatt hour:

$$1 \text{ kWh} = 3.6 \cdot 10^6 \text{ J}$$

Planck constant:

$$h = 6.626 \cdot 10^{-34} \text{ J s}$$

Speed of light in vacuum:

$$c = 2.998 \cdot 10^8 \text{ m s}^{-1}$$

Elementary charge:

$$e = 1.6022 \cdot 10^{-19} \text{ C}$$

Electrical power:

$$P = \Delta E \times I$$

Power efficiency:

$$\eta = P_{\text{obtained}} / P_{\text{applied}}$$

Planck-Einstein relation:

$$E = hc / \lambda$$

Ideal gas equation:

$$pV = nRT$$

Gibbs free energy:


$$G = H - TS$$

$$\Delta_r G^\circ = -RT \ln K^\circ$$

$$\Delta_r G^\circ = -n F E_{\text{cell}}^\circ$$

$$\Delta_r G = \Delta_r G^\circ + RT \ln Q$$

Reaction quotient  $Q$  for a reaction



$$Q = \frac{[\text{C}]^c [\text{D}]^d}{[\text{A}]^a [\text{B}]^b}$$

Henderson–Hasselbalch equation:

$$\text{pH} = \text{p}K_a + \log \frac{[\text{A}^-]}{[\text{AH}]}$$

Nernst–Peterson equation:

$$E = E^\circ - \frac{RT}{zF} \ln Q$$

where  $Q$  is the reaction quotient of the reduction half-reaction

$$\text{at } T = 298 \text{ K}, \frac{RT}{F} \ln 10 \approx 0.059 \text{ V}$$

Beer–Lambert law:

$$A = \epsilon l c$$

Rate laws in integrated form:

- Zero order:

$$[\text{A}] = [\text{A}]_0 - kt$$

- First order:

$$\ln[\text{A}] = \ln[\text{A}]_0 - kt$$

- Second order:

$$1/[\text{A}] = 1/[\text{A}]_0 + kt$$

Half-life for a first order process:

$$t_{1/2} = \ln 2 / k$$

Number average molar mass  $M_n$ :

$$M_n = \frac{\sum_i N_i M_i}{\sum_i N_i}$$

Mass average molar mass  $M_w$ :

$$M_w = \frac{\sum_i N_i M_i^2}{\sum_i N_i M_i}$$

Polydispersity index  $I_p$ :

$$I_p = \frac{M_w}{M_n}$$

## Note

The unit of molar concentration is either “M” or “ $\text{mol L}^{-1}$ ”:

$$1 \text{ M} = 1 \text{ mol L}^{-1}$$

$$1 \text{ mM} = 10^{-3} \text{ mol L}^{-1}$$

$$1 \text{ }\mu\text{M} = 10^{-6} \text{ mol L}^{-1}$$

## Periodic table

|                   |  |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |  |
|-------------------|--|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|--|
| 1                 |  |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   | 18                |                   |                   |                   |  |
| 1<br>H<br>1.008   |  | 2                 |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   | 13                |                   | 14                | 15                | 16                | 17                | 2<br>He<br>4.003  |                   |  |
| 3<br>Li<br>6.94   |  | 4<br>Be<br>9.01   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   | 5<br>B<br>10.81   | 6<br>C<br>12.01   | 7<br>N<br>14.01   | 8<br>O<br>16.00   | 9<br>F<br>19.00   | 10<br>Ne<br>20.18 |  |
| 11<br>Na<br>22.99 |  | 12<br>Mg<br>24.31 |                   | 3                 | 4                 | 5                 | 6                 | 7                 | 8                 | 9                 | 10                | 11                | 12                | 13<br>Al<br>26.98 | 14<br>Si<br>28.09 | 15<br>P<br>30.97  | 16<br>S<br>32.06  | 17<br>Cl<br>35.45 | 18<br>Ar<br>39.95 |                   |  |
| 19<br>K<br>39.10  |  | 20<br>Ca<br>40.08 | 21<br>Sc<br>44.96 | 22<br>Ti<br>47.87 | 23<br>V<br>50.94  | 24<br>Cr<br>52.00 | 25<br>Mn<br>54.94 | 26<br>Fe<br>55.85 | 27<br>Co<br>58.93 | 28<br>Ni<br>58.69 | 29<br>Cu<br>63.55 | 30<br>Zn<br>65.38 | 31<br>Ga<br>69.72 | 32<br>Ge<br>72.63 | 33<br>As<br>74.92 | 34<br>Se<br>78.97 | 35<br>Br<br>79.90 | 36<br>Kr<br>83.80 |                   |                   |  |
| 37<br>Rb<br>85.47 |  | 38<br>Sr<br>87.62 | 39<br>Y<br>88.91  | 40<br>Zr<br>91.22 | 41<br>Nb<br>92.91 | 42<br>Mo<br>95.95 | 43<br>Tc<br>-     | 44<br>Ru<br>101.1 | 45<br>Rh<br>102.9 | 46<br>Pd<br>106.4 | 47<br>Ag<br>107.9 | 48<br>Cd<br>112.4 | 49<br>In<br>114.8 | 50<br>Sn<br>118.7 | 51<br>Sb<br>121.8 | 52<br>Te<br>127.6 | 53<br>I<br>126.9  | 54<br>Xe<br>131.3 |                   |                   |  |
| 55<br>Cs<br>132.9 |  | 56<br>Ba<br>137.3 | 57-71             | 72<br>Hf<br>178.5 | 73<br>Ta<br>180.9 | 74<br>W<br>183.8  | 75<br>Re<br>186.2 | 76<br>Os<br>190.2 | 77<br>Ir<br>192.2 | 78<br>Pt<br>195.1 | 79<br>Au<br>197.0 | 80<br>Hg<br>200.6 | 81<br>Tl<br>204.4 | 82<br>Pb<br>207.2 | 83<br>Bi<br>209.0 | 84<br>Po<br>-     | 85<br>At<br>-     | 86<br>Rn<br>-     |                   |                   |  |
| 87<br>Fr<br>-     |  | 88<br>Ra<br>-     | 89-103            | 104<br>Rf<br>-    | 105<br>Db<br>-    | 106<br>Sg<br>-    | 107<br>Bh<br>-    | 108<br>Hs<br>-    | 109<br>Mt<br>-    | 110<br>Ds<br>-    | 111<br>Rg<br>-    | 112<br>Cn<br>-    | 113<br>Nh<br>-    | 114<br>Fl<br>-    | 115<br>Mc<br>-    | 116<br>Lv<br>-    | 117<br>Ts<br>-    | 118<br>Og<br>-    |                   |                   |  |

|                   |                   |                   |                   |               |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |
|-------------------|-------------------|-------------------|-------------------|---------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|
| 57<br>La<br>138.9 | 58<br>Ce<br>140.1 | 59<br>Pr<br>140.9 | 60<br>Nd<br>144.2 | 61<br>Pm<br>- | 62<br>Sm<br>150.4 | 63<br>Eu<br>152.0 | 64<br>Gd<br>157.3 | 65<br>Tb<br>158.9 | 66<br>Dy<br>162.5 | 67<br>Ho<br>164.9 | 68<br>Er<br>167.3 | 69<br>Tm<br>168.9 | 70<br>Yb<br>173.0 | 71<br>Lu<br>175.0 |
| 89<br>Ac<br>-     | 90<br>Th<br>232.0 | 91<br>Pa<br>231.0 | 92<br>U<br>238.0  | 93<br>Np<br>- | 94<br>Pu<br>-     | 95<br>Am<br>-     | 96<br>Cm<br>-     | 97<br>Bk<br>-     | 98<br>Cf<br>-     | 99<br>Es<br>-     | 100<br>Fm<br>-    | 101<br>Md<br>-    | 102<br>No<br>-    | 103<br>Lr<br>-    |



## Definition of GHS statements

The GHS hazard statements (H-phrases) associated with the materials used are indicated in the problems. Their meanings are as follows.

### Physical hazards

H225 Highly flammable liquid and vapor.  
H226 Flammable liquid and vapor.  
H228 Flammable solid.  
H271 May cause fire or explosion; strong oxidizer.  
H272 May intensify fire; oxidizer.  
H290 May be corrosive to metals.

### Health hazards

H301 Toxic if swallowed.  
H302 Harmful if swallowed.  
H304 May be fatal if swallowed and enters airways.  
H311 Toxic in contact with skin.  
H312 Harmful in contact with skin.  
H314 Causes severe skin burns and eye damage.  
H315 Causes skin irritation.  
H317 May cause an allergic skin reaction.  
H318 Causes serious eye damage.  
H319 Causes serious eye irritation.  
H331 Toxic if inhaled.  
H332 Harmful if inhaled.  
H333 May be harmful if inhaled.  
H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
H335 May cause respiratory irritation.  
H336 May cause drowsiness or dizziness.  
H351 Suspected of causing cancer.  
H361 Suspected of damaging fertility or the unborn child.  
H371 May cause damage to organs.  
H372 Causes damage to organs through prolonged or repeated exposure.  
H373 May cause damage to organs through prolonged or repeated exposure.

### Environmental hazards

H400 Very toxic to aquatic life.  
H402 Harmful to aquatic life.  
H410 Very toxic to aquatic life with long-lasting effects.  
H411 Toxic to aquatic life with long-lasting effects.  
H412 Harmful to aquatic life with long-lasting effects.

## Chemicals

## For all problems

| Chemicals   | Labeled as      | GHS hazard statements |
|---|-----------------|-----------------------|
| Deionized water in:<br>- Wash bottle (bench)<br>- Plastic bottle (bench)<br>- Plastic canister (hood) | Deionized Water | Not hazardous         |
| Ethanol, in a wash bottle   | Ethanol         | H225, H319            |
| Sample of white wine, 300 mL in amber plastic bottle  | Wine sample     | H225, H319            |

## For problem P1

| Chemicals   | Labeled as          | GHS hazard statements                    |
|---|---------------------|--|
| 4-nitrobenzaldehyde, 1.51 g in amber glass vial                                 | 4-nitrobenzaldehyde | H317, H319                               |
| Eluent A, 20 mL in glass vial   | Eluent A            | H225, H290, H304, H314, H319, H336, H410 |
| Eluent B, 20 mL in glass vial   | Eluent B            | H225, H290, H304, H314, H319, H336, H410 |
| Oxone <sup>®</sup> (potassium peroxomonosulfate salt), 7.87 g in plastic bottle | Oxone <sup>®</sup>  | H314                                     |
| Sample of 4-nitrobenzaldehyde for TLC   | TLC standard        | H317, H319                               |

## For problem P2

| Chemicals   | Labeled as                    | GHS hazard statements |
|---|-------------------------------|-----------------------|
| 1 M potassium thiocyanate solution, 20 mL in plastic bottle       | KSCN 1 M                      | H302+H312+H332, H412  |
| 0.00200 M potassium thiocyanate solution, 60 mL in plastic bottle | KSCN 0.00200 M                | Not hazardous         |
| 1 M perchloric acid solution, 10 mL in plastic bottle             | HClO <sub>4</sub>             | H290, H315, H319      |
| 0.00200 M iron(III) solution, 80 mL in plastic bottle             | Fe(III) 0.00200 M             | Not hazardous         |
| 0.000200 M iron(III) solution, 80 mL in plastic bottle            | Fe(III) 0.000200 M            | Not hazardous         |
| 0.3% hydrogen peroxide solution, 3 mL in amber glass bottle       | H <sub>2</sub> O <sub>2</sub> | Not hazardous         |

## For problem P3

| Chemicals   | Labeled as                                      | GHS hazard statements  |
|---|---|------------------------|
| 0.01 M iodine solution, 200 mL in brown plastic bottle                              | <b>I<sub>2</sub></b>                            | H372                   |
| 0.03 M sodium thiosulfate solution, 200 mL in plastic bottle                        | <b>Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub></b> | Not hazardous          |
| 1 M NaOH solution, 55 mL in plastic bottle  | <b>NaOH</b>                                     | H290, H314             |
| 2.5 M sulfuric acid solution, 80 mL in plastic bottle                               | <b>H<sub>2</sub>SO<sub>4</sub></b>              | H290, H315, H319       |
| 0.5 M potassium iodide solution, 25 mL in plastic bottle                            | <b>KI</b>                                       | H372                   |
| Potassium iodate, <i>ca</i> 100 mg (exact mass written on the label), in glass vial | <b>KIO<sub>3</sub></b>                          | H272, H315, H319, H335 |
| Starch solution, 25 mL in plastic bottle  | <b>Starch</b>                                   | Not hazardous          |

**Equipment  
For all problems**

| Personal equipment   | Quantity  |
|--|-----------|
| Pipette filler bulb  | 1         |
| Safety goggles   | 1         |
| 1 L plastic bottle for organic waste, labeled “ <b>Organic waste</b> ” | 1         |
| Paper towels   | 15 sheets |
| Precision wipers   | 30 sheets |
| Spatula (large)  | 1         |
| Spatula (small)  | 1         |
| Stopwatch  | 1         |
| Pencil   | 1         |
| Eraser   | 1         |
| Black pen  | 1         |
| Felt-tip pen for glassware   | 1         |
| Ruler  | 1         |

| Shared equipment              | Quantity  |
|-------------------------------|---|
| UV lamp for TLC visualization | 2 per lab   |
| Colorimeter                   | 5 per lab   |
| Gloves                        | All sizes (S, M, L, XL) available upon request to a lab assistant |
| Ice bucket                    | 1 per lab   |

**For problem P1**

| Personal equipment  | Quantity |
|---|----------|
| Laboratory stand with:  | 1        |
| - Clamp holder with small clamp   | 2        |
| - Clamp holder with large clamp   | 1        |
| Erlenmeyer flask with ground joint, 100 mL  | 1        |
| Erlenmeyer flask with ground joint, 50 mL   | 1        |
| Reflux condenser  | 1        |
| Hotplate stirrer  | 1        |
| Crystallizing dish  | 1        |
| Magnetic stirring bar   | 1        |
| Suction flask   | 1        |
| Büchner funnel with rubber adapter  | 1        |
| Zipped bag with 3 pieces of filter paper  | 1        |
| Petri dish  | 1        |
| TLC elution chamber, labeled “ <b>TLC elution chamber</b> ”                           | 1        |
| Zipped bag with 3 TLC plates (with fluorescence indicator), labeled with Student Code | 1        |
| TLC graduated spotters (in the Petri dish)  | 4        |
| Plastic tweezers  | 1        |
| Glass rod (found in the drawer with glass pipettes)                                   | 1        |
| Graduated cylinder, 25 mL   | 1        |
| Beaker, 150 mL  | 2        |
| Plastic powder funnel (found in the cupboard with glassware)                          | 1        |



|   |   |
|---|---|
| Disposable/Beral plastic pipettes   | 2 |
| Amber glass vial, for TLC sample, 1.5 mL, with stopper, labeled <b>C</b> and <b>R</b> | 2 |
| Pre-weighed amber glass vial, 10 mL, with stopper, labeled with <b>Student Code</b>   | 1 |
| Magnetic stirring bar retrieve  | 1 |

**For problem P2**

| Personal equipment                      | Quantity |
|---|----------|
| Volumetric pipette, 10 mL               | 1        |
| Graduated pipette, 10 mL                | 3        |
| Graduated pipette, 5 mL                 | 3        |
| Test tube stand                         | 1        |
| Test tube                               | 15       |
| Test tube stopper                       | 7        |
| Colorimeter cuvette, path length 1.0 cm | 2        |
| Beaker, 100 mL                          | 2        |
| Disposable/Beral plastic pipettes       | 15       |

**For problem P3**

| Personal equipment                     | Quantity  |
|--|-----------|
| Laboratory stand with burette clamp    | 1         |
| Burette, 25 mL                         | 1         |
| Glass transfer funnel                  | 1         |
| Erlenmeyer flask, 100 mL               | 3         |
| Erlenmeyer flask, 250 mL               | 3         |
| Beaker, 150 mL                         | 1         |
| Beaker, 100 mL                         | 2         |
| Volumetric flask, 100 mL, with stopper | 1         |
| Volumetric pipette, 50 mL              | 1         |
| Volumetric pipette, 25 mL              | 1         |
| Volumetric pipette, 20 mL              | 1         |
| Graduated cylinder, 25 mL              | 1         |
| Graduated cylinder, 10 mL              | 1         |
| Graduated cylinder, 5 mL               | 1         |
| Disposable/ Beral plastic pipettes     | 3         |
| Parafilm                               | 20 sheets |

| Problem P1<br>13% of total | Question | Yield | Purity | TLC | P1.1 | P1.2 | Total |
|----------------------------|----------|-------|--------|-----|------|------|-------|
|                            | Points   | 12    | 12     | 8   | 2    | 3    | 37    |
|                            | Score    |       |        |     |      |      |       |

### Problem P1. Greening the oxidation of nitrobenzaldehyde

For decades, chemists have tried to replace harmful reagents in oxidation processes in order to reduce hazardous waste. In this problem, you will use potassium peroxomonosulfate (Oxone<sup>®</sup>) because it only produces non-toxic and non-polluting sulfate salts. The reaction is performed in a mixture of water and ethanol, which are classified as green solvents.

**Your task is to perform the oxidation of 4-nitrobenzaldehyde, to recrystallize the product, to compare TLC eluents and to check the purity of the product using TLC.**

**Note:** Ethanol waste and eluent must be disposed of in the “Organic waste” bottle.

### Procedure

#### I. Oxidation of 4-nitrobenzaldehyde

- Mix **20 mL of water** and **5 mL of ethanol**.
- Insert the magnetic bar in the 100 mL **ground-joint** Erlenmeyer flask.
- Transfer the pre-weighed **1.51 g of 4-nitrobenzaldehyde** into the Erlenmeyer flask. Add **all** of the **water/ethanol** mixture prepared previously. Clamp the Erlenmeyer flask to the stand. Start stirring the mixture, then add the pre-weighed **7.87 g of Oxone<sup>®</sup>**.
- Attach the reflux condenser by loosening the large clamp and adjusting the ground joints (see Figure 1). Raise your HELP card. A lab assistant will come to turn on the water and set the hotplate.
- Heat the reaction mixture enough to maintain a gentle reflux (*approximately* 1 drop refluxing per second) for **45 minutes**. The mark on the heater corresponds to the setting which should provide a gentle reflux. Monitor your reaction to ensure your heat setting maintains a gentle reflux. Do not change the heat setting unnecessarily.




Figure 1. Setup for heating the reaction mixture under reflux

6. Then turn off the heating on the hotplate stirrer. Remove the Erlenmeyer from the hot plate and let the reaction mixture cool for **10 minutes**. After **10 minutes**, place the Erlenmeyer with the reaction mixture in an ice-water bath using the crystallizing dish. Let the reaction mixture stand for **10 minutes** in the ice-water bath.
7. Set up a vacuum filtration apparatus (see Figure 2) using a Büchner funnel, a filter paper and a suction flask. Ensure that you have secured the set-up to the laboratory stand. Raise your HELP card. A lab assistant will come and assist with setting up the vacuum.




Figure 2. Setup for the vacuum filtration

8. Wet the filter paper with water and ensure that it covers all the holes of the Büchner funnel.
9. Pour the suspension of the crude product into the Büchner funnel and apply vacuum. Wash the solid thoroughly with deionized water (at least **4×20 mL**).
10. Leave the vacuum on for 5 minutes to dry the product. Disconnect the vacuum. Use the small spatula to transfer a small spatula tip amount of **the crude product** to the 1.5 mL amber glass vial labeled **C**. Close the vial and save it for part III (TLC).
11. Transfer **ALL of the remaining solid** into the 50 mL ground-joint Erlenmeyer flask.
12. Discard the filtrate in the “Organic waste” bottle and wash both the suction flask and the Büchner funnel with ethanol and water. Use the “Organic waste” bottle to dispose of ethanol waste.

## II. Recrystallization of the product

- Mix **9 mL of water** and **21 mL of ethanol**.
- Recrystallize the crude product in the 50 mL ground-joint Erlenmeyer flask adding an **appropriate amount** of the **water/ethanol** mixture using the same setup as before for the reflux heating (see Figure 1). Raise your HELP card. A lab assistant will come to turn on the water and set the hotplate. If needed, add more **water/ethanol** solvent from the top of the condenser.
- Once the product has crystallized, use the same cooling and drying procedure as described previously (I.7 to I.10) to collect the solid. Use the small spatula to transfer a small spatula tip of amount the **recrystallized product** into the 1.5 mL amber glass vial, labeled **R**. Close the vial and save it for part III (TLC).

4. **Carefully** transfer the **ALL** of the rest of the recrystallized solid in the pre-weighed vial labeled with your **Student Code**. Close the vial.
5. Discard the filtrate in the “Organic waste” bottle and raise your HELP card. A lab assistant will come to turn off the water of the condenser.

### III. TLC analysis

1. Prepare the elution chamber. Load the elution chamber with *approximately* 0.5 cm depth of **eluent A**. Cover the elution chamber with a Petri dish. Wait for the eluent to saturate the atmosphere in the elution chamber.
2. Prepare your samples. You have a sample of **4-nitrobenzaldehyde** in an amber glass vial labeled **TLC standard** (referred as **S** on the TLC). You have a small sample of your crude product (vial **C**) and your recrystallized product (vial **R**) in two other amber glass vials. Add *approximately* **1 mL of ethanol** to each of the vials to dissolve the samples.
3. Prepare your TLC plate. Use a pencil to draw the start line carefully 1 cm above the bottom of the plate. Mark the positions where you will spot the 3 samples. Label them **S** (TLC standard), **C** (Crude product) and **R** (Recrystallized product), as depicted in Figure 3. On the top left of the plate, write your **Student Code**. On the top right of the plate, write the eluent you use (**Eluent A** on the first plate, then **Eluent B** when you repeat at step 7). Spot the three samples on the plate, using capillary spotters.




Figure 3. TLC plate preparation

4. Perform the TLC analysis. Using tweezers, insert the TLC plate into the elution chamber and cover it with the Petri dish. Let the eluent reach approximately 1 cm below the top of the plate. Using tweezers, remove the plate, mark the eluent front with a pencil and let the plate air-dry.
5. Place the TLC plate under one of the UV lamps on the common bench. With a pencil, carefully circle all the visible spots.
6. Discard the eluent into the “Organic waste” bottle.
7. Repeat steps 1, 3, 4, 5, and 6 with **eluent B**.
8. Place your plates in the zipped bag with your Student Code.

Results of your TLC analysis (complete the schemes with your results). You may use these drawings to represent the results of your TLC plates, which may help you answer the following questions. The drawn representation of your results (below) will not be graded.

| Eluent A   |  |  | Eluent B   |  |  |
|--|--|--|--|--|--|
| <div style="position: absolute; top: 0; left: 0; right: 0; border-top: 1px dashed black; height: 10px;"></div> <div style="position: absolute; bottom: 0; left: 0; right: 0; border-bottom: 1px dashed black; height: 10px;"></div> <div style="position: absolute; bottom: 0; left: 10%; width: 20%; text-align: center;">S</div> <div style="position: absolute; bottom: 0; left: 35%; width: 20%; text-align: center;">C</div> <div style="position: absolute; bottom: 0; left: 60%; width: 20%; text-align: center;">R</div> |  |  | <div style="position: absolute; top: 0; left: 0; right: 0; border-top: 1px dashed black; height: 10px;"></div> <div style="position: absolute; bottom: 0; left: 0; right: 0; border-bottom: 1px dashed black; height: 10px;"></div> <div style="position: absolute; bottom: 0; left: 10%; width: 20%; text-align: center;">S</div> <div style="position: absolute; bottom: 0; left: 35%; width: 20%; text-align: center;">C</div> <div style="position: absolute; bottom: 0; left: 60%; width: 20%; text-align: center;">R</div> |  |  |

At the end of the examination, your lab supervisor will pick up the following items:

- Glass vial labeled with your **Student Code** containing your recrystallized product;
- TLC plates A and B in zipped bag labeled with your **Student Code**.

|  |                          |  |
|--|--------------------------|--|
| Submitted items  |                          |  |
| <b>Recrystallized product</b>  | <input type="checkbox"/> |  |
| <b>TLC plate A</b>   | <input type="checkbox"/> |  |
| <b>TLC plate B</b>   | <input type="checkbox"/> |  |
| <b>Signatures</b> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="text-align: center; width: 45%;"> <hr style="border: 0; border-top: 1px solid black; margin-bottom: 5px;"/>             Student           </div> <div style="text-align: center; width: 45%;"> <hr style="border: 0; border-top: 1px solid black; margin-bottom: 5px;"/>             Lab Supervisor           </div> </div> |                          |  |

**Questions**

1. Propose a structure for the final organic product from the reaction of 4-nitrobenzaldehyde and Oxone<sup>®</sup>.

2. Based on your results on the TLC analysis, answer the following questions.

- Which eluent is better to follow the reaction progress?

☐ **A** ☐ **B**

- The crude product (C) contains traces of 4-nitrobenzaldehyde.

☐ **True** ☐ **False**

- The recrystallized product (R) contains traces of 4-nitrobenzaldehyde.

☐ **True** ☐ **False**

| Problem<br>P2<br>14% of<br>total | Question | Calibration | Iron<br>determination | P2.1 | P2.2 | P2.3 | Stoichiometry<br>determination | P2.4 | P2.5 | Total |
|----------------------------------|----------|-------------|-----------------------|------|------|------|--------------------------------|------|------|-------|
|                                  | Points   | 10          | 6                     | 3    | 4    | 3    | 9                              | 3    | 2    | 40    |
|                                  | Score    |             |                       |      |      |      |                                |      |      |       |

### Problem P2. The iron age of wine

Iron is an element which can naturally be found in wine. When its concentration exceeds 10 to 15 mg per liter, iron(II) oxidation into iron(III) may lead to quality loss, through the formation of precipitates. It is therefore necessary to assess the iron content of the wine during its production.


Given the very low concentration of iron species, a colored complex of iron(III) with thiocyanate  $\text{SCN}^-$  as a ligand is used to quantify the iron amount, through spectrophotometric measurements.

**Your task is to determine the total iron concentration of the white wine provided, using spectrophotometry, and to determine the stoichiometry of the thiocyanate – iron(III) complex.**


### WARNING

- In this task, you are provided with two iron(III) solutions and two potassium thiocyanate solutions of different concentrations. Be very careful not to confuse them.
- Once the solutions are ready for spectrophotometric measurements, record the absorbance no later than one hour after the addition of thiocyanate.
- When you need a colorimeter, raise your HELP card. A lab assistant will give you a specific colorimeter (labelled) for all of your measurements. Be sure to record the code of the colorimeter you use in your results table where “colorimeter code” is indicated. You will have the exclusive use of this colorimeter for up to 15 minutes. The lab assistant will take the colorimeter back as soon as you have finished or when the 15 minutes are over. If no colorimeters are available at the precise moment you need them, you will be added to a waiting-list. Try to use the same colorimeter throughout your experiment but if the same colorimeter is not available when you are ready to measure absorbance, go with the available colorimeter and be sure to record the code.
- Instructions for the colorimeter on the next page.
- You can use the colorimeter **ONLY three times** (i.e. for only three 15 min intervals) for this problem.

## Instructions for the use of the colorimeter



- Plug in the colorimeter.
- Check that “Absorbance” is highlighted. If not, turn the selection wheel until a dashed line appears around “Absorbance” and then press the OK button.
- Turn the selection wheel until a dashed line appears around the desired wavelength (470 nm). Press the OK button.
- Place the cuvette with *approximately* 3 cm (height) of the blank solution in the cuvette. Be careful to choose the correct orientation for the cuvette (look at the orientation scheme on the colorimeter, the beam is in the direction of the yellow arrow, see figure below), and to push the cuvette down fully. Close the lid.
- Turn the selection wheel until a dashed line appears around “Absorbance” and then press the OK button. Using the selection wheel, highlight “Calibration” and press the OK button.
- Wait until the display reads 0.00 (or –0.00).
- Place the cuvette with *approximately* 3 cm (height) of the analyzed solution in the cuvette. Close the lid.
- Read the absorbance value.





## I. Determination of the iron content in the wine

In this part, you will need the 0.000200 M iron(III) solution and the 1 M potassium thiocyanate solution.

### Procedure

1. Prepare 6 test tubes by adding the volumes of solutions, as described in the table below.


| Tube #                        | 1      | 2      | 3      | 4      | 5       | 6       |
|-------------------------------|--------|--------|--------|--------|---------|---------|
| 0.000200 M iron(III) solution | 1.0 mL | 2.0 mL | 4.0 mL | 6.0 mL |         |         |
| 1 M perchloric acid solution  | 1.0 mL | 1.0 mL | 1.0 mL | 1.0 mL | 1.0 mL  | 1.0 mL  |
| Wine                          |        |        |        |        | 10.0 mL | 10.0 mL |
| Hydrogen peroxide solution    |        |        |        |        | 0.5 mL  | 0.5 mL  |
| Deionized water               | 9.5 mL | 8.5 mL | 6.5 mL | 4.5 mL |         | 1.0 mL  |

2. Stopper and shake the test tubes.
3. Add **1.0 mL of 1 M potassium thiocyanate** solution in tubes **1, 2, 3, 4 and 5 only**. Do **not** add potassium thiocyanate to tube **6**. Stopper and shake all solutions.
4. When all the tubes are prepared, raise your HELP card to get a colorimeter from a lab assistant.
5. Prepare the colorimeter using the procedure described previously (see page 16). Set the wavelength at 470 nm. Use deionized water for the blank.
6. Record the absorbance of each tube (**1 to 6**) at 470 nm. Report the results in the table below. Raise your HELP card to return the colorimeter.

| Tube #  | 1  | 2  | 3  | 4  | 5 | 6 |
|---|----|----|----|----|---|---|
| Absorbance (at 470 nm)  |    |    |    |    |   |   |
| Analytical concentration of $\text{Fe}^{3+}$ in the tube<br>$c(\text{Fe}^{3+}) / \mu\text{M}$ | 16 | 32 | 64 | 96 |   |   |
| Colorimeter code  |    |    |    |    |   |   |

## Questions

- Plot the absorbance  $A$  of tubes **1** to **4** as a function of the analytical concentration of  $\text{Fe}^{3+}$  in the tube.



- To show which values you selected as relevant for determining your calibration curve, mark the relevant values with an “X” in the table below. Leave the box in the table blank if the value is an outlier.

| Tube #   | 1 | 2 | 3 | 4 |
|--|---|---|---|---|
| Absorbance values used in the calibration curve marked with an “X” |   |   |   |   |

2. Using the calibration plot and the data you have chosen, **draw** the calibration curve (it should be a straight line) on the graph. Show your calculations and determine the analytical concentration (in  $\mu\text{mol L}^{-1}$ ) of  $\text{Fe}^{3+}$  in tube **5** using the calibration curve.

$$c(\text{Fe}^{3+})_{\text{TUBE 5}} = \underline{\hspace{2cm}} \mu\text{mol L}^{-1}$$

*If you could not calculate  $c(\text{Fe}^{3+})$ , the value  $c(\text{Fe}^{3+}) = 50 \mu\text{mol L}^{-1}$  can be used in the rest of the problem.*

3. Calculate the mass concentration, in mg per liter, of iron in the studied white wine.

$$c_{\text{m}}(\text{iron}) = \underline{\hspace{2cm}} \text{mg L}^{-1}$$

## II. Determination of the complex stoichiometry

In this part, you will need the 0.00200 M iron(III) solution and the 0.00200 M potassium thiocyanate solution.

### Procedure

In part I of this problem, we used the color of the iron(III)-thiocyanate complex to determine the concentration of iron in the sample of wine. In part II of this problem, you will determine the stoichiometry of the  $[\text{Fe}_a(\text{SCN})_b]^{(3a-b)+}$  complex (coordination of water is not shown), where  $a$  and  $b$  are integers no greater than 3.

You are provided with the following aqueous solutions for this part:

- 0.00200 M iron(III) solution (already acidified) (80 mL)
- 0.00200 M potassium thiocyanate solution (80 mL)

You also have test tubes (with stoppers that you can wash and dry), graduated pipettes, a spectrophotometer cuvette, a colorimeter (upon request), and any other labware on your bench that you think is useful.

1. Record the volumes and the mole fraction you will use to determine the stoichiometry of the complex by spectrophotometric measurements. *You don't have to fill all the columns but the more points you have, the more precise your curve might be.* Consider time and precision when you prepare your tubes. Calculate the molar fraction of iron(III) in each tube, using the following formula.

$$x(\text{Fe}^{3+}) = \frac{V_{\text{Fe(III)}}}{V_{\text{Fe(III)}} + V_{\text{SCN}^-}}$$


| Tube #   | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 |
|--|---|---|---|----|----|----|----|----|----|
| Volume of 0.00200 M iron(III) solution<br>$V_{\text{Fe(III)}} / \text{mL}$           |   |   |   |    |    |    |    |    |    |
| Volume of 0.00200 M potassium thiocyanate solution<br>$V_{\text{SCN}^-} / \text{mL}$ |   |   |   |    |    |    |    |    |    |
| Molar fraction in iron(III) $x(\text{Fe}^{3+})$                                      |   |   |   |    |    |    |    |    |    |
| Absorbance (at 470 nm)   |   |   |   |    |    |    |    |    |    |
| Colorimeter code   |   |   |   |    |    |    |    |    |    |

2. Prepare the tubes you will use for your absorbance measurements. When all the tubes you want to use are ready, raise your HELP card to get a colorimeter from a lab assistant.
3. Prepare the colorimeter using the procedure described previously (see page 16). Set the wavelength at 470 nm. Use deionized water for the blank.

4. Record the absorbance of each tube at 470 nm. Report the results in the table above (on the previous page).

### Questions

4. Plot the absorbance  $A$  of the tubes as a function of the molar fraction of iron(III)  $x(\text{Fe}^{3+})$ .



5. Based on the results of the experiments you carried out, determine the stoichiometry of the complex  $[(\text{Fe})_a(\text{SCN})_b]^{(3a-b)+}$ .


$a =$  \_\_\_\_\_

$b =$  \_\_\_\_\_

| Problem<br>P3<br>13% of<br>total | Question | Titration<br>I | Titration<br>II | Titration<br>III | P3.1 | P3.2 | P3.3 | P3.4 | P3.5 | Total |
|----------------------------------|----------|----------------|-----------------|------------------|------|------|------|------|------|-------|
|                                  | Points   | 10             | 10              | 8                | 4    | 4    | 2    | 2    | 2    | 42    |
|                                  | Score    |                |                 |                  |      |      |      |      |      |       |

### Problem P3. Wine for keeping

Sulfur dioxide,  $\text{SO}_2$ , is used as a preservative in wine. When  $\text{SO}_2$  is added to wine, it can react with water leading to bisulfite ions,  $\text{HSO}_3^-$ , and protons,  $\text{H}^+$ . Bisulfite can also be converted to sulfite,  $\text{SO}_3^{2-}$ , by the loss of a second proton.



These three different forms of sulfur dioxide in water can react with chemicals in wine such as acetaldehyde, pigments, sugars, etc. forming products P. The total concentration of sulfur dioxide is the sum of the concentration of the “free” forms ( $\text{SO}_2$ ,  $\text{HSO}_3^-$  and  $\text{SO}_3^{2-}$ ) and P.

The preservative concentration is regulated because sulfites and sulfur dioxide can be harmful to some people. In the EU, the maximum total sulfur dioxide content is set at  $100 \text{ mg L}^{-1}$  for red wine and  $150 \text{ mg L}^{-1}$  for white or rosé.

**Your task is to determine the total sulfur dioxide concentration of white wine by iodometric titration.**

### Procedure

#### I. Standardization of the sodium thiosulfate solution

1. You are given a sample of *approximately 100 mg of pure potassium iodate*  $\text{KIO}_3$ . The exact mass is written on the label of the vial. Report it in the table below.
2. Prepare **100.00 mL of potassium iodate** solution in the 100 mL volumetric flask, using the **whole sample of solid potassium iodate and deionized water**. This solution is called S.
3. In a 100 mL Erlenmeyer flask, add:
  - **20.00 mL of solution S** with a volumetric pipette;
  - **5 mL of the potassium iodide solution (0.5 M)**, with a 5 mL graduated cylinder;
  - **10 mL of the sulfuric acid solution (2.5 M)** with a 10 mL graduated cylinder.
4. Swirl the Erlenmeyer flask, cover it with Parafilm and keep it in the cupboard for at least five minutes.
5. Fill the burette with the thiosulfate solution using a beaker. Titrate the contents of the Erlenmeyer flask with constant swirling. When the liquid turns pale yellow, **add ten drops of the starch solution** and keep titrating until the solution becomes colorless. Record the titration volume  $V_1$ .
6. Repeat the procedure (steps 3-5) as needed.

| Mass of potassium iodate<br>(report the value on the label) |            |
|---|------------|
|   |            |
| Analysis n°   | $V_1$ / mL |
| 1   |            |
| 2   |            |
| 3   |            |
|   |            |
|   |            |
| <b>Reported value <math>V_1</math> / mL</b>                 |            |

## II. Standardization of the iodine solution

1. With a volumetric pipette, transfer **25.00 mL of the iodine solution labeled  $I_2$**  into a 100 mL Erlenmeyer flask.
2. Titrate the content of the Erlenmeyer flask with the sodium thiosulfate solution. When the liquid turns pale yellow, **add ten drops of the starch** solution and keep titrating until the solution becomes colorless. Record the titration volume  $V_2$ .
3. Repeat the procedure (steps 1-2) as needed.

| Analysis n°                                 | $V_2$ / mL |
|---|------------|
| 1   |            |
| 2   |            |
| 3   |            |
|   |            |
|   |            |
| <b>Reported value <math>V_2</math> / mL</b> |            |

**III. Determination of total sulfur dioxide**

1. With a volumetric pipette, transfer **50.00 mL of wine** into a 250 mL Erlenmeyer flask.
2. Add **12 mL of the sodium hydroxide solution (1 M)**, with a 25 mL graduated cylinder. Cover the flask with Parafilm, swirl the content then let it stand for at least **20 minutes**.
3. Add **5 mL of the sulfuric acid solution (2.5 M)**, and *approximately 2 mL of starch solution* using a graduated disposable plastic pipette.
4. Titrate the content of the Erlenmeyer flask with the iodine solution in the burette, until a dark color appears and persists for at least 15 seconds. Record the titration volume  $V_3$ .
5. Repeat the procedure (steps 1-4) as needed.

| Analysis n°                                 | $V_3$ / mL |
|---|------------|
| 1   |            |
| 2   |            |
| 3   |            |
|   |            |
|   |            |
| <b>Reported value <math>V_3</math> / mL</b> |            |



**Questions**

1. Write down the balanced equations of all the reactions occurring during the standardization of the sodium thiosulfate solution.

2. Calculate the molar concentration of the sodium thiosulfate solution. The molar mass of potassium iodate is  $M(\text{KIO}_3) = 214.0 \text{ g mol}^{-1}$ .

$$c(\text{S}_2\text{O}_3^{2-}) = \underline{\hspace{2cm}} \text{ mol L}^{-1}$$

*If you could not calculate  $c(\text{S}_2\text{O}_3^{2-})$ , the value  $c(\text{S}_2\text{O}_3^{2-}) = 0.0500 \text{ mol L}^{-1}$  can be used in the rest of the problem.*

3. Calculate the molar concentration of the iodine solution.

$$c(\text{I}_2) = \underline{\hspace{2cm}} \text{mol L}^{-1}$$

*If you could not calculate  $c(\text{I}_2)$ , the value  $c(\text{I}_2) = 0.00700 \text{ mol L}^{-1}$  can be used in the rest of the problem.*

4. Write down the equation of the reaction between iodine  $\text{I}_2$  and sulfur dioxide  $\text{SO}_2$ , assuming that sulfur dioxide is oxidized into sulfate ions  $\text{SO}_4^{2-}$ .

5. Calculate the mass concentration, in mg per liter, of total sulfur dioxide in the wine. The molar mass of sulfur dioxide is  $M(\text{SO}_2) = 64.1 \text{ g mol}^{-1}$ .

$$c_{\text{m}}(\text{SO}_2) = \underline{\hspace{2cm}} \text{mg L}^{-1}$$

## PENALTIES

| Incident #     | Student signature | Lab supervisor signature |
|----------------|-------------------|--------------------------|
| 1 (no penalty) |                   |                          |
| 2              |                   |                          |
| 3              |                   |                          |
| 4              |                   |                          |
| 5              |                   |                          |